



## Sim.STREAM

### Top-quality telemedicine

In today's operating theaters, telemedicine plays an increasingly important role. OT designers are faced with new challenges: structure distances need to be bridged, so as to enable seamless, interactive collaboration between the theater, offices, training facilities or external institutions.

With Sim.STREAM, SIMEON offers you the perfect solution for meeting these requirements. Sim.STREAM gives you the capability to transmit a video signal over the network simultaneously to multiple PCs in superior HD quality and at a resolution of 1920 x 1080 pixels. At the same time, the video signal and single still frames can both be easily recorded to a USB flash drive.



PUSHING TECHNOLOGY TO EXCELLENCE



- Simultaneous retrieval of the video stream from different sources
- Individual hierarchy and password assignment
- Simultaneous control of Sim.STREAM by front keypad and PC software
- Up to 6 video and 4 audio inputs

Sim.STREAM	
Video inputs	2x HD-SDI, 2x DVI, 2x Y/C or CVBS
Video outputs	1x HD-SDI, 1x DVI
Audio inputs	4x stereo
Simultaneous record and stream	yes
USB port	1x (USB 3.0-compatible)

Tolerance ±10%; technical specifications are subject to change



S.I.M.E.O.N. Medical GmbH & Co. KG  
 In Grubenäcker 18  
 78532 Tuttlingen, Germany  
 Phone: +49 7461 90068-0  
 Fax: +49 7461 90068-900  
[www.simeonmedical.com](http://www.simeonmedical.com)  
[info@simeonmedical.com](mailto:info@simeonmedical.com)



PUSHING TECHNOLOGY TO EXCELLENCE