



Sim.CAM 4K
More Details. Better Treatment.



PUSHING TECHNOLOGY TO EXCELLENCE

The new dimension in video images

Fascinating images – with clear, sharp contours and strikingly improved color. Details that Full HD was unable to depict are now clear and distinct with 4K ULTRA HD.



The Sim.CAM 4K operating panel



External camera

Perfectly integrated, highly flexible, easy-to-clean surfaces




4K technology

Improves the sharpness, color, and contrast of your surgical images

Easy to operate

Camera functions adjustable right on the camera arm



Use 4K in surgical settings for detail-optimized imagery. Its exceptional picture quality provides finely structured images with more accurate details, and thus more brilliant visual representations.

Sim.CAM 4K	
Resolution	8,847,360 Pixel
4K ULTRA HD	4096 x 2160
Zoom	12x optical 12x digital
Image reorientation	> 360°
Video outputs	HDMI 1.4 6G SDI

Tolerance ±10 %; technical specifications are subject to change



PUSHING TECHNOLOGY TO EXCELLENCE