



The new Dimension in Video Images

Wireless




NEW: Sim.CAM 4K Wireless* –

More Details. Better Treatment.

Fascinating images—with clear, sharp contours and strikingly improved color. Details that Full HD was unable to depict are now clearly and distinctly visualized in 4K ULTRA HD. If needed, a cable version of Sim.CAM 4K as a separate arm solution is available.

Sim.CAM 4K and Sim.CAM 4K Wireless*	
Resolution	8,847,360 pixels
4K Ultra HD	4096 x 2160
Zoom	12x optical 12x digital
Image positioning	> 360°, mechanically and electronically
Video outputs	HDMI 1.4, DVI

Technical specifications are subject to change
* Available starting 3rd quarter/2018



Use 4K ULTRA HD in surgical settings for detail-optimized imagery. The high image quality offers you finely structured images with amazing details and thus a brilliant visual representation.

NEW: Sim.CAM HD Wireless –

More Details. Great Flexibility.

With the new Sim.CAM HD Wireless you obtain a premium HD camera that sends images to your monitor without additional wiring, thanks to the new wireless technology. The integrated HD camera transfers all data sharply and precisely from every selected position. If needed, it is also available as a cable-based solution.

Sim.CAM HD and Sim.CAM HD Wireless	
Resolution	2,000,000 pixels
Full HD	1,920 x 1,080
Aspect ratio	16:9
Zoom	10x optical 12x digital
Image positioning	> 360°, mechanically and electronically
Video outputs	HDMI, DVI

Technical specifications are subject to change